

Paleobotany, Micropaleontology, & Mineralogy Collection Sam Noble Oklahoma Museum of Natural History

Scanning Supplies for Procedures


The following list is a general list of types of supplies in use or available to be used (it is not intended to be a complete detailed account of all equipment); if you need further information please consult the Collection Manager.

- Computer with Keyboard & Mouse
- Epson Perfection 1650 scanner



(Photo of Epson Perfection 1650 Scanner)

- Photoshop (used to access scanner through Epson TWAIN driver and save images)
- Epson TWAIN driver (used to control scanner through Photoshop)
- SNOMNH Slide Template Transparency



(Photo of SNOMNH Slide Template Transparency)

[Printed on overhead transparency from file: SNOMNH_Slide_Scan_Temp.pdf]

- Alcohol (used to periodically clean scanner glass)
- Kimwipes EXL (used to periodically clean scanner glass and SNOMNH Slide Template Transparency)



(Photo of Kimwipes)

